



LUSEAUX®

HYDRION QT TEST PAPER

MEASURES QUATERNARY AMMONIUM CONCENTRATIONS

GENERAL DESCRIPTION: A test paper in a convenient dispenser which gives instant readings of the concentration of quaternary ammonium compounds in sanitizing solutions. Match points on the color chart indicate 100-200-300-400 ppm of "quat" concentration. Will allow the user to detect exhaustion of solutions that should be replaced and will reduce the use of excessive and wasteful amounts of "quat" compound. Each dispenser comes with a color chart and enough paper for well over 120 tests.

APPLICATION: Open dispenser and remove paper roll from foil wrapping. Place the roll back into the dispenser and replace dispenser top. Tear approximately 2" of Hydrion QT test paper from the dispenser and dip into the solution to be tested for **10 SECONDS EXACTLY**. Make sure to stay away from any foam and do **NOT** stir the test strip as it will cause a high reading. Remove after the 10 seconds are up and **immediately** match to the color chart. Refer to sanitizer label for information on required concentrations.

Always consult your local health department for more information regarding sanitizing regulations for your area.

LUSEAUX® Laboratories, Inc.

DETERGENTS • CLEANERS • SANITIZERS

Gardena, CA 90247 • 800-266-1555

Kingman, AZ 86401 • 800-895-4583